

Effective January 1, 2007 GE Healthcare (Amersham Biosciences) is no longer an authorized dealer of KODAK Molecular Imaging Systems Autoradiography and Electron Micrography film products.

Kodak Molecular Imaging Systems Autoradiography film products including:

BioMax MR film	(Maximum Resolution)
BioMax MS film	(Maximum Sensitivity)
BioMax Light film	(Maximum Performance for Chemiluminescence)
BioMax XAR film	(Excellent General Purpose film)
X-Omat Blue (XB) film	(Economical film)
BioMax MS intensifying screens	(Maximum Sensitivity in conventional screen format)
BioMax cassettes	(Maximize image Resolution & Clarity)
BioMax TranScreen intensifying screens	(Maximum Results with an Innovative Screen)
Darkroom Illuminator LED Safelight	(Maximize your image quality)

(To install the Kodak Illuminator safelight simply screw the unit into the existing light socket in your darkroom, this safelight can be used with all Autoradiography and EM films.)

are all available through

PerkinElmer Life and Analytical Sciences
www.perkinelmer.com

KODAK BioMax Performance Guide

Label	Maximum Sensitivity	High Sensitivity High Resolution	Maximum Resolution
High Energy Radioisotopes (i.e. ^{32}P & ^{125}I)	BioMax TranScreen HE screen/BioMax MS film	BioMax TranScreen LE screen/BioMax MS film	BioMax MR film
Low and medium energy radioisotopes (i.e. ^3H , ^{14}C , ^{35}S , ^{33}P & ^{59}Fe)	BioMax TranScreen LE screen/BioMax MS film		BioMax MR film
Fluorography	BioMax TranScreen LE screen/BioMax MS film*		
Chemiluminescence	BioMax Light film		BioMax Light film

* Radioactive particle must be accessible to the screen