

# KODAK Darkroom Illuminator LED Safelight

**NEW!**

The New KODAK Darkroom Illuminator is designed for better autoradiography and electron micrography images



## Improve Your Images

- Produces better quality images - no fogging from worn or damaged filters
- Better visibility in your darkroom
- Cluster of 20 light emitting diodes in the film safe red spectrum

## No Maintenance

- Self-contained unit ready for immediate use
- Easy to install into a standard light socket
- Eliminates need to replace bulbs and filters

## Affordable

- Costs less than traditional darkroom lamps, bulbs, and filters
- Full three year warranty

## Cat #

## Description

113 2752 KODAK Darkroom Illuminator, North American socket E26, 115V

897 4917 KODAK Darkroom Illuminator, European socket E27, 230V

155 6208 KODAK Darkroom Illuminator, UK socket BA-22, 230V

For more information or technical assistance visit our website at [www.kodak.com/go/biomax](http://www.kodak.com/go/biomax) or email to [sis-support@kodak.com](mailto:sis-support@kodak.com) or call +1-203-786- 5657

**KODAK SCIENTIFIC IMAGING SYSTEMS**

